



Glooko[®] Uploader

Instructions for use – for Glooko users

Contents

Find your device	4-5
General information	6
Product description.....	6
Intended use.....	6
Supported devices.....	6
Warnings.....	6
Contraindications	6
Instructions for use	7
Patients	8
Sign into your MyGlooko dashboard at https://my.glooko.com	8
Download the software installation file	8
Locate the installation file and run the file.....	8
Healthcare providers.....	9
Request a link to download the Glooko Uploader software.....	9
Download the software installation file	9
Locate the installation file and run the file.....	9
Patients and healthcare providers	10
Open the Glooko Uploader software.....	10
Sign into your Glooko account.....	10
Features in Glooko Uploader	11
Upload diabetes data	12
View diabetes data.....	13
Compatible devices and transmission methods	14

Find your device



Glucose meter



CGM



Insulin pump

Abbott

	FreeStyle Freedom Lite.....	14
	FreeStyle InsuLinx	14
	FreeStyle Lite.....	14
	FreeStyle Optium Neo	14
	FreeStyle Precision Neo.....	14
	Precision Xtra.....	14
	Precision Xceed.....	14

	FreeStyle Libre.....	14
	FreeStyle Libre 2.....	14

Acon Laboratories

	On Call® Express.....	14
	On Call® Express II.....	14
	On Call® Express Voice	14
	On Call® Sure.....	14
	On Call® Sure Sync	14

AgaMatrix

	WaveSense Jazz™	14
--	-----------------------	----

Animas

	OneTouch Ping (meter).....	15
--	----------------------------	----

Animas

	OneTouch Ping (pump).....	20
	Vibe	20

Arkray

	GlucoCard Shine Connex	15
	GlucoCard Shine Express	15
	GlucoCard Shine XL	15
	ReliOn Premier Blu.....	15
	ReliOn Premier Classic.....	15
	ReliOn Premier Voice.....	15

Ascensia (Bayer)

	Breeze 2.....	15
	Contour	15
	Contour Care	15
	Contour Fit.....	15
	Contour Link	15
	Contour Next.....	15
	Contour Next (Bluetooth)	15
	Contour Next One	15
	Contour Plus.....	15
	Contour Plus Blue	15
	Contour Plus Elite.....	15
	Contour Plus One.....	15
	Contour next EZ	15
	Contour next Link.....	15
	Contour next Link 2.4	15
	Contour next USB.....	15
	Contour TS.....	15
	Contour XT	15
	Contour USB.....	15

Beurer

	GL 44.....	16
	GL 48.....	16
	GL 49.....	16

Bionime

	Rightest GM700SB.....	16
--	-----------------------	----

Biosys

	Aria	16
--	------------	----

Dexcom

	G4 Platinum	19
	G4 Platinum with Share	19
	G5	19
	Touchscreen Receiver (G5/G6).....	19

Ekoweb

	GlucoX (TD-4183).....	16
	GlucoX Pro (TD-4183).....	16

Embrace

	TALK.....	16
--	-----------	----

Equil/Wellion

	MICRO-pump.....	21
--	-----------------	----

ForaCare

	Diamond PRIMA	16
--	---------------------	----

GlucoRx

	Nexus.....	16
	Nexus Mini	16
	Q.....	16

GlucoX

	TD-4183.....	16
--	--------------	----

i-SENS

	alphacheck professional	16
	CareSens BGM/TEE2.....	16
	CareSens Dual.....	16
	CareSens N	16
	CareSens N POP	16
	CareSens N Premier	16
	CareSens N Voice.....	16
	No Coding Plus	16
	TEE2+	16

Infopia	
 Finetest Lite	16

Insulet	
 Omnipod®	21
 Omnipod® Dash™ System	21

Intuity	
 POGO® Automatic™	17

LifeScan	
 OneTouch Select Plus	17
OneTouch Ultra.....	17
OneTouch Ultra Plus Reflect.....	17
OneTouch Ultra2	17
OneTouch UltraEasy.....	17
OneTouch UltraMini	17
OneTouch Verio	17
OneTouch Verio Flex.....	17
OneTouch Verio IQ.....	17
OneTouch Verio Reflect.....	17
OneTouch Verio Sync.....	17

Medtronic	
 630G.....	21
 640G.....	21
670G.....	21

Menarini	
 GlucoFix Tech	17
GlucoMen areo.....	17
GlucoMen areo 2K.....	17
GlucoMen LX2.....	17

Nipro	
 4SURE Smart	17
4SURE Smart Duo.....	17

Prodigy	
 AutoCode.....	18

Roche	
 Accu-Chek Active.....	18
Accu-Chek Aviva Connect.....	18
Accu-Chek Aviva Combo.....	19
Accu-Chek Aviva Expert.....	19
Accu-Chek Aviva Nano.....	19
Accu-Chek Aviva Plus Black	19
Accu-Chek Aviva Plus Silver	19
Accu-Chek Compact Plus.....	19
Accu-Chek Guide	18
Accu-Chek Guide Me	18
Accu-Chek Instant.....	18
Accu-Chek Mobile (USB)	18
Accu-Chek Nano.....	19
Accu-Chek Performa	19
Accu-Chek Performa Nano	19

Roche	
 Accu-Chek Aviva Insight	22
Accu-Chek Performa Insight.....	22
Accu-Chek Aviva Solo.....	21
Accu-Chek Guide Solo	21
Accu-Chek Spirit Combo	22

Sanofi	
 BGStar	18

Taidoc	
 TD-4141	18
TD-4183	18
TD-4255	18
TD-4257(v1)	18
TD-4257(v2)	18
TD-4280	18
TD-4285	18
TD-4289	18

Tandem	
 t:flex.....	21
t:slim.....	21
t:slim G4.....	21
t:slim X2	21

Terumo	
 Medisafe Fit Smile	20
 Medisafe WITH.....	20

ViCentra	
 Kaleido.....	21

Trividia (Nipro)	
 TRUE METRIX	17
TRUE METRIX AIR.....	17
TRUEresult.....	17
TRUEyou mini.....	17

Wellion	
 Calla Mini.....	18

XPÉR Technology	
 A1 (TD-4183).....	18

Ypsomed	
 Mylife Aveo.....	18

General information

Product description

Glooko Uploader® is indicated for use by individuals or healthcare professionals in the home or healthcare facilities to transfer predefined data from home monitoring devices to a server database. The device can be used as a stand-alone product or in combination with the Glooko System.

Intended use

The Glooko System is a diabetes management platform that can be used in-clinic or remotely by Health Care Professionals and patients. It allows viewing of patient data from compatible medical devices and other health monitoring devices. It also allows Health Care Professionals to support the treatment of diabetes through analysis and monitoring of a patient's diabetes therapy and Care Program.

Supported devices

To get started, individuals will need:

1. A compatible blood glucose meter, insulin pump or continuous glucose monitor (CGM);
2. A PC or Mac computer with one of the following operating systems (OS):
 - Windows 10 (32-bit and 64-bit)
 - Windows 11 (32-bit and 64-bit)
 - macOS 11.0: Big Sur
 - macOS 12.0: Monterey
 - macOS 13.0: Ventura
3. the Glooko Uploader software;
4. A Glooko account; and
5. The syncing or device download cable supplied with the supported diabetes device.

Warning!

Glooko does not measure, interpret, or make decisions on the data it conveys nor is it intended to provide automated treatment decisions or be used as a substitute for professional judgment. All medical diagnosis and treatment are to be performed under the supervision and oversight of an appropriate healthcare provider.

If you have concerns about your current diabetes diagnosis and treatment, please contact your healthcare provider.

In the event where the Glooko System is malfunctioning, please discontinue use of the product and contact Glooko Support.

Contraindications


None.

Instructions for use – Glooko Uploader

You will need to create a Glooko account to use Glooko Uploader. The data uploaded using Glooko Uploader will synchronise to your Glooko account, which can be accessed using your Glooko login information. Internet (WiFi or mobile data network) access is necessary for your data to synchronise to the Glooko server.

Patients:

- Learn how to install and upload data using Glooko Uploader on the next page.



Healthcare providers:

- Learn how to install and upload data using the Glooko Uploader on page 9.



Patients


In order to access the MyGlooko dashboard and download the Glooko Uploader software, you will need to have a Glooko account. You can sign up for Glooko using the Glooko mobile app for iOS or Android. You can also contact **Glooko Support** (<https://support.glooko.com>) for assistance in creating an account. Once you have created a Glooko account, please follow the instructions below.

1 Sign into your MyGlooko dashboard at <https://my.glooko.com>

- Enter your email address and password.
- Tap Sign In.

GET HELP SIGNING INTO YOUR GLOOKO ACCOUNT

- **FORGOT PASSWORD:**
If you are unable to remember your password to sign in, tap **Forgot Password** to receive a password reset email. Follow the instructions on the screen. You can also contact Glooko Support for assistance.
- **UNLOCK YOUR ACCOUNT:**
If you enter your password incorrectly five times, your account will be locked. An email with unlock instructions will automatically be sent to the email address registered to your account. However, if you do not receive an email, click **Didn't receive unlock instructions** – you will be redirected to the Unlock page. Follow the instructions on the screen. You can also contact Glooko Support for assistance.



2 Download the software installation file

After logging into your MyGlooko dashboard, you will need to download the installation file.

- Tap **Settings**.
- Scroll to the bottom of the page and click **Get Glooko Uploader**.
- You will be taken to a screen with buttons to install Glooko Uploader on a Mac or PC. Click either the Mac or PC button to download the installation file.



3 Locate the installation file and run the file

Depending on your web browser and operating system, you may need to search for the downloaded installation file or you may see a popup window asking you to run the file.

- If you are using a **PC**, look in your recent downloads folder
- If you are using a **Mac**, look in the downloads folder in Finder

Once you have located the file, double-click on it and follow the installation instructions on the screen.
NOTE: You will need to restart your computer at the end of the installation process.

Healthcare providers

In order to use the Glooko Uploader software, you will need a Glooko provider account and a subscription that supports Glooko Uploader. You will also need a specific link to download Glooko Uploader. Please contact Glooko Support (<https://support.glooko.com>) for assistance.

1 Request a link to download the Glooko Uploader software

Healthcare providers must request a link to download the Glooko Uploader software.

- Please contact Glooko Support (<https://support.glooko.com>) to request a link to download the Glooko Uploader software.

2 Download the software installation file

- When you receive the link to download the software from Glooko, click the link in your email. You will be taken to a screen with buttons to install Glooko Uploader on a Mac or PC. Click either the Mac or PC button to download the installation file.



3 Locate the installation file and run the file

Depending on your web browser and operating system, you may need to search for the downloaded installation file or you may see a popup window asking you to run the file.

- If you are using a **PC**, look in your recent downloads folder
- If you are using a **Mac**, look in the downloads folder in Finder

Once you have located the file, double-click on it and follow the installation instructions on the screen.

NOTE: You will need to restart your computer at the end of the installation process.

Patients and healthcare providers

After installing Glooko Uploader and restarting your computer, locate the Glooko Uploader icon on either your desktop (PC) or in your Applications folder (Mac).

4 Open the Glooko Uploader software

- **Locate** the Glooko Uploader icon on your desktop (PC) or in your Applications folder (Mac).
- **Double-click** on the Glooko Uploader icon to start the program.

NOTE: Make sure you have restarted your computer before starting Glooko Uploader for the first time.



5 Sign into your Glooko account

Once the Glooko Uploader program opens, you will need to sign into your Glooko account.

- Enter your **email or username**
- Enter your **password**
- Tap **Sign In**.

GET HELP SIGNING INTO YOUR GLOOKO ACCOUNT

If you cannot remember your password to sign in, click **Forgot Password** to receive a password reset email. You can also contact Glooko Support for assistance.

Patients and healthcare providers


6 Features in Glooko Uploader

LOG OUT OF YOUR ACCOUNT

- You can log out of your Glooko account by tapping the Profile button  in the bottom right-hand corner of Glooko Uploader.
- Tap **Log out**


VIEW YOUR PROFILE INFORMATION

- You can view your profile information by tapping the Profile button  in the bottom right-hand corner of Glooko Uploader.
- This popup includes:
 - o Your **name**
 - o Your **primary authorised site** (healthcare provider) and the associated **ProConnect code**
 - o Your **email address**
- Tap **Log Out**.



SEE ADDITIONAL INFORMATION ABOUT GLOOKO UPLOADER

- You can view additional information by tapping the Info button  in the bottom left-hand corner of Glooko Uploader.
- This popup includes product support information, product description, and the intended use statement.



Patients and healthcare providers


7 Upload diabetes data

Once you have signed into Glooko Uploader, you can upload your diabetes device data.

CONNECT THE CABLE FOR YOUR DIABETES DEVICE TO THE COMPUTER

After logging in, Glooko Uploader will start searching for a device.

- Connect the appropriate USB cable or infrared (IR) cable for the diabetes device to the computer.
- Once the cable is recognised, wait a few seconds for the cable drivers to be recognised before connecting your device.




CONNECT THE DIABETES DEVICE TO THE CABLE

- After connecting the cable to the computer, connect the diabetes device to the cable.

NOTE: Make sure you use the cable supplied with your diabetes device. If you need a new cable, please contact the manufacturer of your diabetes device.

- After connecting the diabetes device to your computer using the cable supplied with the device, your device will begin uploading data using Glooko Uploader.
- If an error occurs, please follow the instructions on the screen to resolve the error. If you continue to experience difficulties, please contact Glooko Support for assistance.



Patients and healthcare providers

8 View diabetes data


After uploading your diabetes device data using Glooko Uploader, you can view the data at <https://my.glooko.com> by logging into your Glooko account.

- **Patients:**
To learn more about the MyGlooko web dashboard, please refer to [Glooko for personal use – Quick start guide](#) at <https://support.glooko.com>
- **Healthcare providers:**
To learn more about Population Tracker, please refer to: [Glooko for clinics – Quick start guide](#) at <https://support.glooko.com>



Compatible devices and transmission methods

Below are instructions on how to upload diabetes data to Glooko using Glooko® Uploader. Unless explicit upload instructions are provided for the device, it will upload automatically when you connect it to the PC/Mac using its specific USB cable.



Please note that some of the diabetes devices may not be available in your market.

Blood glucose meters – connected via USB cable

Abbott

FreeStyle Freedom Lite

FreeStyle Lite

Abbott FreeStyle USB
(2,5 mm)



FreeStyle InsuLinx
FreeStyle Precision Neo

FreeStyle Optium Neo

Micro USB



Precision Xceed

Precision Xtra

Abbott Xceed
USB



FreeStyle Libre

FreeStyle Libre 2

Micro USB



Compatibility is restricted to
Sweden and Norway by Abbott.

Acon Laboratories

On Call® Express
On Call® Express II
On Call® Express Voice

On Call® Sure
On Call® Sure Sync

On Call USB
(2,5 mm)



Agamatrix

WaveSense Jazz

Agamatrix USB



Blood glucose meters – connected via USB cable

Animas

OneTouch Ping (meter)

Mini USB



1. Connect to your computer.
2. Press the power button on the meter.

Arkray

Glucocard Shine XL

i-Sens USB
(2.5 mm)



Connect and press any button.

Glucocard Shine Connex
Glucocard Shine Express

Micro USB



ReliOn Premier Blu (BLE)
ReliOn Premier Classic
ReliOn Premier Voice

i-Sens USB
(2.5 mm)



Connect and press any button.

Ascensia (Bayer)

Breeze 2

Contour Plus



1. Connect to your computer.
2. Press the power button on the meter.

Contour
Contour link
Contour next EZ

Contour TS
Contour XT

1. Connect to your computer.
2. Press the **M** button on the meter.

Contour next link
Contour next link 2.4

Contour next USB
Contour USB

USB port on meter




Contour Care
Contour Fit
Contour Next
Contour Next (Bluetooth)

Contour Next one
Contour Plus one
Contour Plus Blue
Contour Plus Elite


Micro USB




Blood glucose meters – connected via USB cable



Blood glucose meters – connected via USB cable



Blood glucose meters – connected via USB cable



Blood glucose meters – with infrared connection



Roche (IR)

Accu-Chek Aviva Nano
Accu-Chek Aviva Plus Black
Accu-Chek Aviva Plus Silver
Accu-Chek Nano
Accu-Chek Performa
Accu-Chek Perfoma Nano



ACTiSYS Li, Li2 or Accu-Chek 360°

1. Connect the IR cable to your computer.
2. With the meter off, activate IR transmission by pushing and holding the two arrow buttons until two arrows appear in the display.

Accu-Chek Compact Plus



ACTiSYS Li, Li2 or Accu-Chek 360°

1. Connect the IR cable to your computer.
2. Activate IR transmission by pushing and holding the **S** and **M** buttons below the display until two arrows appear in the display.

Accu-Chek Aviva Combo
Accu-Chek Aviva Expert



ACTiSYS Li, Li2 or Accu-Chek 360°

1. Connect the IR cable to your computer.
2. In the main menu, go to **My data** and select **Data transfer**.

CGMs – connected via USB cable

Dexcom

G4 Platinum
G4 Platinum with Share
G5
Touchscreen Receiver (G5/G6)

Micro USB



If the upload fails to start, disconnect the receiver and select Shutdown in the menu to turn the receiver off. Turn the receiver on again and connect it to the computer.

Blood glucose meters and Insulin pumps - with NFC connection

Terumo
 Medisafe Fit Smile



1. Connect the NFC plate provided by Terumo to your computer.
2. Place the meter on top of the NFC plate.

Please note that the transmission symbol on the back of the remote control needs to be aligned with the NFC symbol on the NFC Plate.

Medisafe WITH



1. Connect the NFC plate provided by Terumo to your computer.
2. Make sure the main menu is showing on the remote control for the pump.
3. Place the remote control on top of the NFC plate.

Please note that the transmission symbol on the back of the remote control needs to be aligned with the NFC symbol on the NFC Plate.

Insulin pumps – with infrared connection

Animas (IR)
 OneTouch Ping Vibe



1. Suspend the pump. (MAIN MENU, scroll to **Suspnd/Resum** and press **OK**; **Suspend** is highlighted).
2. Press **OK**.
3. Place the pump with its back to the Actisys IR wireless download cable.
4. When the transmission is finished: Resume pump (MAIN MENU, scroll to **Suspnd/Resum** and press **OK**; **Resume** is highlighted). Press **OK**.

ACTiSYS Li2 or Li2
 Note! If you are using Windows 8 or 10, you need to use the ACTiSYS Li2 cable.

Compatible infrared cables


ACTiSYS – Li cable  Compatible operating systems: macOS X: 10.14–12.0 Windows: N/A	ACTiSYS – Li2 cable  Compatible operating systems: macOS X: 10.14–12.0 Windows: 8, 10	Accu-Chek – 360° cable (Realtyme)  Compatible operating systems: macOS X: 10.14–12.0 Windows: 8, 10	Roche Accu-Chek Smart Pix 2 cable  Compatible operating systems: macOS X: 10.14–12.0 Windows: 8, 10
---	--	---	--



Need a cable?

Please contact the manufacturer of your specific diabetes device.

Insulin pumps – connected via USB cable



Insulin pumps – connected via infrared cables

Compatible infrared cables


Roche Accu-Chek
Smart Pix cable



Roche Accu-Chek
Smart Pix 2 cable



Roche
Accu-Chek Aviva Insight
Accu-Chek Performa Insight



Smart Pix 2

Note! Connect the handset of your pump with a micro USB cable to the Smart Pix 2.

1. Connect the Smart Pix 2 cable to your computer and wait for the blue indicator on Smart Pix 2 to start flashing slowly.
2. Connect the handset to Smart Pix 2 using a micro USB cable.
3. In the handset menu select **Connect to PC** and press **OK**.
4. Ensure that the pump is within range of the handset.
5. The blue indicator on Smart Pix 2 is on while data is being transmitted from the device to Smart Pix 2.
6. The blue indicator turns itself off and both Smart Pix 2 and the device emit a sound to signal that data is being transmitted from Smart Pix 2 to Glooko Uploader. Do not disconnect Smart Pix 2 yet!
7. Wait for Glooko Uploader to confirm the upload was successful.

Roche
Accu-Chek Spirit Combo



Smart Pix or Smart Pix 2

Note! Smart Pix software version 3.02 or higher is required when using Smart Pix.

1. Connect the Smart Pix to your computer and wait for the blue indicator on the Smart Pix to start flashing slowly.
2. Put the pump in **Data Transfer** mode (stop the pump and enter the **Data Transfer** menu).
3. Place the pump with its IR eye towards the Smart Pix. The blue indicator signals the data is being transmitted to the Smart Pix.
4. Wait for the blue indicator to start flashing slowly again or stop flashing altogether.
5. The data is now being transmitted from the Smart Pix to Uploader. Do not disconnect the Smart Pix yet!
6. Wait for Uploader to signal the upload was successful.



Glooko Inc
579 University Avenue
Palo Alto, CA 94301
USA



Glooko AB
Nellickevägen 20
SE-412 63 Göteborg
SWEDEN

Support
help@glooko.com
<https://support.glooko.com>