

How to Sync Your Pivot™ Pump with Glooko®

Our clinic uses Glooko® to securely review glucose trends from your diabetes devices. Connecting your accounts is the easiest way for us to see your full diabetes picture in one place. You can get connected in just a few minutes using the Pivot™ Diabetes Mobile App.

CONNECTING VIA THE PIVOT DIABETES MOBILE APP

If you already have a Glooko account


Step 1: Install the Pivot™ Diabetes Mobile App

Download the **Pivot™ Diabetes** Mobile App from the Apple App Store or Google Play Store. You can scan the QR codes below with your phone camera to find the app.



Step 2: Register and Connect to your Glooko Account

Once you install the **Pivot™ Diabetes** app on your phone, select **Sign into Your Glooko Account**



Step 3: Log in to Glooko

Use your existing Glooko username (email) and password to log into your Glooko account.

Log in to Glooko


Email

Password

[Show](#)


[Forgot password?](#)

Step 4: Agree to the Terms and Conditions



If you don't have an existing Glooko account

1. Once you install the **Pivot™ Diabetes** app on your phone, select **Sign Up with Glooko**
2. Fill out the registration form
3. Agree to the terms to continue



Tell us about yourself

If you are creating an account on behalf of someone else, please enter their details below.

First Name

Last Name

Date of Birth


Country of Residence

Your country of residence is required for regulatory and data security reasons.

Phone Number (Optional)

If you opt in, your phone number is used to send you information about this ecosystem and your device via SMS.

Already have an account? [Log in](#)



Once you have signed in, continue on to pair your pump to the app and set a basal schedule.

Need Technical Help?

If you need further assistance, please reach out to Modular Medical's Customer Care team via email at customercare@modular-medical.com or via phone at 866-710-1200.